



Impressions: When to use which material?

	Honigum Pro-Light	Honigum Pro-Mono	Honigum Pro-Heavy	Honigum Pro-Putty	Honigum-Light	Honigum-Mono	Honigum-Heavy	Honigum-Putty	StatusBlue
Crown and bridge impressions	■	■	■	■	■	■	■	■	
Inlay and onlay impressions	■	■	■	■	■	■	■	■	
Implant impressions			■				■		
Fixation impressions			■				■		
Functional impressions		■				■			
Pressure-spot indication		■				■			
Situation impression									■
Impressions for the fabrication of temporaries									■
Impressions for the fabrication of vacuum-formed templates									■
Opposing-jaw models									■
Fabrication of models									■